

4Way Glide System

The 4Way Glide is a system particularly suitable for handling and transferring clients in the bed. It is a unique system as it secures controlled transfer in all directions, with the lowest possible use of energy. The system ensures a safe handling of the client, is very user-friendly and simple to place in bed, even while the client is in the bed. The 4Way Glide system can remain in the bed continuously as it has a low impact on pressure-relieving mattresses.

The 4Way Glide system consists of a nylon sheet that is attached to the bed mattress, either with or without different non-slip locking systems, and a glide mattress placed on top. The glide mattress has a low-friction surface on the underside that facilitates movement to all directions and several handles that ensure good grip, better working position and reduced strain. The upper side of the glide mattress is available either in green color cotton or grey color polyurethane (PU) surface.

The nylon sheet is available in several sizes for adaptation to bed size, while the glide mattress is available in both several sizes and configurations, depending on client, application, and caregiver requirements.

4Way Glide LPL

The 4Way Glide LPL is the basic configuration consisting of a nylon sheet with non-slip lock and glide mattress with handles, all adapted to a single size bed 90 x 200 cm. It is used for movement in all directions, positioning in bed and transferring in and out of bed. Turning can be performed manually or with supporting floor or ceiling lift. Maximum user weight: 200 kg




IM140/200LPL	4Way Glide mattress, cotton/polyester, Green	140 x 200 cm
IM140/200PLG	4Way Glide mattress, polyurethane, Grey	140 x 200 cm
IM85/200PLNS	Nylon sheet with non-slip locking system and non-slip along sides	85 x 200 cm
IM85/200LPL	Nylon sheet with non-slip locking system	85 x 200 cm
IM85/2003LPL	Nylon sheet with 3 x non-slip locking system	85 x 200 cm



4Way Glide Twin

The 4Way Glide Twin has a 2-part glide mattress, especially adapted to situations when the client needs to be turned into the side position for a longer period, such as for lower and upper body hygiene and washing, wound care, etc. It is also used in situations when it is difficult to hold the client on the side position, due to pain, dementia, heavy weight, etc. The 4Way Glide Twin can be combined with a floor or ceiling lift as support for turning the client. Maximum user weight: 200 kg.



IM140/200TWIN	2-parted glide mattress, cotton/polyester, Green	140 x 200 cm
IM140/200PLG	2-parted glide mattress, polyurethane, Grey	140 x 200 cm
IM85/200PLNS	Nylon sheet with non-slip locking system and non-slip along sides	85 x 200 cm
IM85/200LPL	Nylon sheet with non-slip locking system	85 x 200 cm
IM85/2003LPL	Nylon sheet with 3 x non-slip locking system	85 x 200 cm



4Way Glide Maxi

The 4Way Glide Maxi system is particularly suitable for handling bariatric clients in the bed. The unique system secures controlled transfer in all direction by using the least amount of strength. The 4way Glide Maxi alternative has a net enforcement built into the glide mattress, allowing higher user weight and wider configuration for larger beds. Maximum user weight: 300 kg



IM180/200LPL	4Way Glide mattress, cotton/polyester, Green	180 x 200 cm
IM110/200LPLNS	Nylon sheet with non-slip locking system and non-slip along sides	110 x 200 cm



4Way Glide Midi

4Way Glide Midi is suitable for clients that are partially self-supporting and can support the movement with their legs. The system is suitable for transferring and positioning in bed. It consists of a short turning mattress that allows movement in all directions and with a locking system and short nylon sheet with non-slip material along the sides. The non-slip locking system allows the client to sit securely on the edge of the bed. Maximum user weight: 200 kg

IM140/150	4Way Glide mattress, cotton/polyester, Green	140 x 150 cm
IM140/150G	4Way Glide mattress, polyurethane, Grey	140 x 150 cm
IM140/150NS	Nylon sheet with non-slip along sides	140 x 150 cm



4Way Glide Mini

The 4Way Glide Mini is a short turning mattress suitable for positioning and transferring self-supporting clients in bed, as well as for turning out of the bed. The nylon sheet is supplied without the non-slip locking system, and allows low-friction movement to all directions. Maximum user weight: 200 kg

IM140/100	4Way Glide mattress, cotton/polyester, Green	140 x 100 cm
IM140/100G	4Way Glide mattress, polyurethane, Grey	140 x 100 cm
IM140/100NS	Nylon sheet with non-slip along sides	140 x 100 cm



Satin Sheet

The Satin Sheet is a series of bed sheets that reduces friction, suitable for clients that have difficulties moving in bed. The Satin Sheet products consist of a cotton sheet with low-friction (satin) material in the center that reduces friction and therefore allows easier movement in bed. Reducing friction in the center of the bed increases the possibility for a self-sufficient client to move more often and thereby reduces the risk of decubitus. Satin Sheet also supports caregiver-assisted positioning. The stitching between cotton and satin material is lowered reducing the risk of decubitus. The Satin Sheet series has numerous variants in terms of sizes and combinations.

IM4113	Satin Sheet, 2-direction	200 x 140 cm
IM4115	Satin Sheet, 4-direction	200 x 140 cm
IM4118	Satin Sheet, draw sheet	200 x 140 cm

