

## Support Belt

The Support Belt is used when the client is being supported to and from a standing position from a wheelchair, chair, armchair, bed, etc. It is also used for walking support and exercise. The specially designed locking system allows the application of the belt without twisting the back and the non-slip inner surface prevents the belt from sliding up. The handles allow for caregiver support and the diagonal position provides a more natural and ergonomic grip.



IM8010509	Support belt max 150 kg, 4 handles Small	9 x 105 cm
IM8010513	Support belt max 150 kg, 5 handles Small	13 x 105 cm
IM8013013	Support belt max 150 kg, 6 handles Medium	13 x 130 cm
IM8015013	Support belt max 150 kg, 8 handles Large	13 x 150 cm
IM8017013	Support belt max 150 kg, 11 handles X-Large	13 x 170 cm



## Bed String

Bed String with handles is suitable for the self-sufficient client for sitting up in bed, transferring in and out of bed, from bed to wheelchair, turning to side position, dressing, etc. Bed String is attached to the foot end of the bed and the client uses the handles to move into position. The loop allows attachment to most types of bed frames. Maximum user weight: 150 kg



IM429	Bed String with 8 handles	7.5 x 170 cm
-------	---------------------------	--------------



## Multiglide Glove

The Multiglide Glove is suitable for small adjustments of positioning and for smaller transfers, for instance adjustments of shoulder and hip position. It is also very suitable for control of pressure when using pressure relief cushions or mattresses. The material on the external surface minimizes the friction for easy insertion of the hands underneath the patient. The compact size and light weight makes it easy to carry along.



IM35	Multiglide Gloves, pair	20 x 50 cm
------	-------------------------	------------



## SmartMove

SmartMove supports transferring between sitting and standing position, as well as moving in and out of bed. SmartMove allows the caregiver a secure grip and provides a better working position and reduced strain during transfers. The non-slip inner surface prevents the pad from sliding out of position while the low-friction outer surface provides easy application and sliding function. Maximum user weight: 150 kg



IM425	SmartMove 2 handles	21 x 60 cm
-------	---------------------	------------