

Transfer Mattress and Cushion

The Transfer Mattress is developed for regular turning or mobilizing of bed bound clients, and is also suitable for horizontal transfer situations, for instance between beds, operation and x-ray tables, etc. The Transfer Cushion is used for turning and positioning in bed as well as transferring the client to a chair or wheelchair. These are constructed as a continuous padded nylon loop with cotton outer surface and nylon inner surface for smooth and low-friction glide function. It is easy to insert and remove from under the client.



IM60/50	Transfer Cushion	60 x 50 cm
IM60/90	Transfer Cushion	60 x 90 cm
IM90/90	Transfer Cushion	90 x 90 cm
IM70/180	Transfer Mattress	70 x 180 cm
IM70/195	Transfer Mattress	70 x 195 cm
IM60/195	Transfer Mattress	60 x 195 cm



HipTurner

HipTurner is used for turning the client to side position and creates space for placing a bed pan, washing and hygiene, wound care and dressing. HipTurner is easy to apply and remove after use, and provides a better working position and reduced strain for the caregiver. Turning can be performed manually or with support of floor or ceiling lift. Maximum user weight: 350 kg



IM462	HipTurner	Cushion: 40 x 4.5 x 20 cm Glide sheet: 25 x 80 cm
-------	-----------	--



WiseMover

WiseMover is a sliding system that is fixed to the bed and can remain placed under the client. It is suitable as support for positioning and transferring of bedridden clients that are exposed to frequent x-ray situations. The sliding function facilitates positioning and placing the x-ray cassette under the client.

IM441	WiseMover sliding system	85 x 200 cm
-------	--------------------------	-------------

