

Model FG 95
ropimex®-Mobile screen



ropimex®-Mobile screen

The mobile screen '95 is ready for use in a few seconds and quickly back in park position. Its ease of motion, minimum weight, and high stability guarantee simplified use. The mobile screen equipped with one or two telescopic arms is ideal for mobile screening, e.g. for patient accommodation in corridors, in multi-bed rooms, in the emergency department or on the gynaecology seat.



Mobile screen with folded-down telescopic arms.



Mobile screen in hospital.

Mobile screen '95 with two telescopic arms.



FG 95



r

FG 95

ropimex®-Trolley

FG 95




FG 95 Fahrgestell

165x60x40 cm (height/depth/width) [66x24x16"], 4 rollers with a diameter of 100 mm [4"], of which 2 rollers can be locked diagonal in position. Aluminium anodised, glass-fibre reinforced polyamide, corrosion-proof, disinfectable.

LA Electrically conductive rollers available at extra cost.

UKB Universal clasp for attachment of a second telescopic arm. All telescopic arms of the model series RTS and RTI can be inserted in the trolley FG 95.

UKB



Recommended telescopic arms for FG 95:
RTI 2.2 and RTI 1.6

RTI 2.2




RTI 1.6



We recommend the following equipment:

With one telescopic arm:
FG 95 trolley, RTI 2.2 telescopic arm, swivelling and folding down, extended 220 cm [88"], curtain TCS 212 (230x140 cm), [92x56"]

With two telescopic arms:
FG 95 trolley, RTI 2.2 telescopic arm, swivelling and folding down, extended 220 cm [88"], UKB, universal clasp for use of a second telescopic arm, RTI 1.6, swivelling and folding down, extended 160 cm [64"], curtains:
TCS 212 (230x140 cm), [92x56"]
TCS 152 (170x140 cm), [68x56"]



Mobile screen '95 with two telescopic arms.

In park position.

r

r