



ancona lie low

The lie-low bed for optimal prophylaxis from falling out of the bed.



ancona lie low, with a height adjustment of 21–61 cm lying surface with ripoplan



ancona lie low, with a height adjustment of 21–81 cm lying surface with ripoplan

The extremely low lying-height of only 21 cm of the ancona lie-low bed nearly excludes, especially in conjunction with the roll-off mattress, injuries due to falling out of the bed. Therefore even active patients are protected, as they will not climb out of the bed and possibly harm themselves by that.

The central locking for the rollers in a 2+2-system, the large 100 mm steering rollers, the integrated wooden siderails and the lifter space have proved their value in everyday care and can also be found in the ancona lie-low bed. The patient's rehabilitation is supported by the electrically adjustable lying surface with lowered head/leg positions. The security of the bed is ensured by several locks.

On request, the ancona lie-low bed is also available with an adjustable range from 21 to 81 cm and this version fulfils all normative demands. The height adjustment up to 81 cm ensures ideal and back friendly working conditions for care personnel.


Specifications – ancona lie low

Total weight:	145 kg
External dimension:	105 x 222 cm
Safe capacity::	220 kg
Maximum person weight::	185 kg
Height adjustment:	21 – 61 cm oder 21– 81 cm
Lifter space:	> 15 cm
Side rail height:	38 cm with wooden slats, 38 cm with ripoplan
Dimension of lying surface:	90 x 200 cm
Special sizes:	100 cm width, 180 – 220 cm length

Motors


Water-resistant to IPX4, high-grade anti-kink , pull relief and anti-kink protection, fuse in the mains plug and emergency lowering device.

Subject to technical modification.



4-section
Electrical adjustment of the back and leg part with automatic triple function.

Electrical height adjustment



On demand with lowered head/leg position (Trendelenburg).

Extras:

You can find descriptions of the accessories on page 46 – 47.

Roll-down mattress for lie low bed	○
Lifting pole with triangle grip	○
Grab rail	○
Raised leg position 90°	○
Raised leg position stretched	○
Bed lamp model D	○
Bed lamp model E	○
Bed extension for the length of 220 cm	○
Wooden slats with degree of tension	●
Deposit for the bed clothes	○
Infusion holder for lifting pole	○
Infusion holder double (fixed)	○
Infusion holder quadruple (telescopic)	○
Integrated wooden side rails	●
Mattresses	○
Mattress compensation	○
Mattress extension for a length of 220 cm	○
Mechanical bed version	○
Integrated head rest for 5-section lying surface	○
Low entrance height	○
Bumpers for side rail attachment	○
Bumpers for side rails	○
ripolux®- lying surface	○
rioplan- lying surface	○
Covering panel of the lying surface frame	●
Side rail attachment	○
Side rails in special height *	○
Long splitted side rails	○
Short splitted side rails	○
Special width/length	○
Control box and hand control for Trendelenburg	●
Individual disabling function	●
Special mattresses	○
Steel lying surface	○
Transport and storage system	○
Universal clamp with modular accessories	○
Lower leg part with manual adjustment (raster)	○
Wall protection roll (vertical/horizontal)	○
Wall protection roll (gummed)	○
Wall protection buffer head part/ lateral	○
Central braking	●

● Standard ○ Extras

* not possible with splitted side rails

Safety guarantees

- Mains isolation system
- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Disabling function
- Hand control with single fault safety
- Special spiral mains cable
- Pull relief at the mains cable with anti-kink protection
- Water-resistant to IPX4
- Distortion-free, highly stable side rails

