

livorno

The versatile model

Livorno is the quick-change artist among the object beds, because nearly no other bed can be adjusted so perfectly to the individual taste and to external general conditions.

The livorno adapts to any ambience: whether to already existing furnishings or to new living spaces under construction.

The wide range of design and colour variations as well as the technical furnishing possibilities of the livorno object bed offer optimal conditions for more individuality for the furnishing.

Livorno – the convincing all-rounder!







Side rail model
long splitted



livorno

The versatile model

The livorno healthcare bed is extremely comfortable. The fully automatic lying surface offers a quadfunction option. To suit individual needs, the lying surface can also be equipped with the innovative ripolux® system. livorno takes care of healthy, comfortable lying in any position.

A range of furniture for any ambience. No other healthcare furniture line from Bock comes in as many different designs and with as many technology options.

Specifications – livorno


Total weight:	109 kg
External dimension:	105 x 212 cm
Safe capacity::	200 kg
Maximum person weight::	165 kg
Height adjustment:	37 – 83 cm
Lifter space:	15 cm
Side rail height:	39 cm with wooden slats 34 cm with ripolux®
Dimension of lying surface:	90 x 200 cm
Special sizes:	80 – 140 cm width, 180 – 220 cm length
Special equipment:	threeparted, continous side rails (picture on the right side) distance piece on the head end, also available as a telescopic model (without picture) lateral distance piece (picture below right side)

Motors	Water-resistant to IPX4, high-grade anti-kink , pull relief and anti-kink protection, fuse in the mains plug and emergency lowering device.
---------------	---

Variations of equipment

End panels	Type I / II / III / V / VI / VII / VIII
Side rails	continous wooden side rails splitted wooden side rails

Subject to technical modification.



Extras:

You can find the descriptions of the accessories on page 46 – 47.

Roll-down mattress for lie low bed	
Lifting pole with triangle grip	○
Grab rail	○
Raised leg position 90°	○
Raised leg position stretched	○
Bed lamp model D	○
Bed lamp model E	○
Bed extension for the length of 220 cm	○
Wooden slats with degree of tension	●
Deposit for the bed clothes	○
Infusion holder for lifting pole	○
Infusion holder double (fixed)	○
Infusion holder quadruple (telescopic)	○
Integrated wooden side rails	●
Mattresses	○
Mattress compensation	○
Mattress extension for a length of 220 cm	○
Mechanical bed version	○
Integrated head rest for 5-section lying surface	○
Low entrance height	○
Bumpers for side rail attachment	○
Bumpers for side rails	○
ripolux®- lying surface	○
ripoplan- lying surface	○
Covering panel of the lying surface frame	●
Side rail attachment	○
Side rails in special height *	○
Long splitted side rails	○
Short splitted side rails	○
Special width/length	○
Control box and hand control for Trendelenburg	
Individual disabling function	○
Special mattresses	○
Steel lying surface	○
Transport and storage system	○
Universal clamp with modular accessories	○
Lower leg part with manual adjustment (raster)	○
Wall protection roll (vertical/horizontal)	○
Wall protection roll (gummed)	○
Wall protection buffer head part/ lateral	○
Central braking	○

● Standard ○ Extras

* not possible with splitted side rails

Safety guarantees

- Mains isolation system
- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Disabling function
- Hand control with single fault safety
- Special spiral mains cable
- Pull relief at the mains cable with anti-kink protection
- Water-resistant to IPX4
- Distortion-free, highly stable side rails

