

2 Move Transfer Board

The 2 Move transfer board is specially designed for horizontal transfer of patients between, for instance, a bed, operating table, X-ray table, bath or washing table, etc. Using the 2 Move transfer board, friction is minimized, differences in height can be overcome and strain on the personnel is reduced. The 2 Move is easy to transport and store and the long model can be folded for less space. Using the adapted Multiglide cover increases sliding function as well as external protection and hygiene control. Adapted disposable tube foil is available for further requirements on anti-bacterial control.



IM8050090	2 Move transfer board max 150 kg	50 x 90 cm
IM5050180	2 Move transfer board max 150 kg Foldable at 90 cm	50 x 180 cm
IM8150180	2 Move transfer board max 150 kg Non-foldable	50 x 180 cm
IM4044	PU-Nylon waterproof cover for 2 Move	50 x 90 cm
IM4045	PU-Nylon waterproof cover for 2 Move	50 x 180 cm
IM4030	Disposable tube foil cover for 2 Move Roll with 75 pcs	60 x 200 cm
IM4031	Disposable tube foil cover for 2 Move Roll with 150 pcs	60 x 100 cm



Mini Board

The Mini Board is used for turning and transferring between bed and operating table, x-ray table, etc as well as for transfer situations from the floor. Maximum user weight: 150 kg



IM8012454	Mini Board 7 handles	54 x 124 cm
IM8020054	Mini Board 11 handles	54 x 200 cm

Scanner Arm Support

Scanner Arm Support cushion provides a comfortable arm positioning for clients that have difficulties maintaining the arms fixed in position during CT and MR scanning. The cushion provides a more comfortable position and prevents movement, increasing the quality of the scanning. The surface can easily be wiped off with disinfectant solution.



IM500	Scanner Arm Support	30 x 25 x 52 cm
-------	---------------------	-----------------

